

Pressure Mapping



Pressure mapping is a specialised measurement technology used to measure and visualize the contact pressure distribution between the human body and a supporting surface and equipment interface, e.g. person, chair or sling. Care & Independence commission independent pressure mapping experts to conduct such trials to ascertain sling performance and help identify areas of risk. The subsequent scientific data insight has enabled Care & Independence to develop solutions and vastly improve upon the areas which indicate tissue viability risks, pain or other health concerns to the equipment user.

115


RESULTS FOR GLOVE™ IN-CHAIR DIVIDED LEG (GIDL)

The data from these tests is impressive. After ten minutes, the performance of the GLOVE™ In-Chair Divided Leg sling with the Polyair® cushion is 11 mmHg lower than without the GIDL in-situ. This shows that the addition of the GLOVE™ In-Chair Divided Leg can actually be better for a user.

In the baseline mapping tests where the subject was sat clothed direct upon both cushion types, a clear low pressure result was returned as indicated by the expanse of blue colouring. [Fig. 1 & 2]

In the two minute seating test, the **GIDL/Polyair® combination shows a decrease to peak pressure.** With the Integrity® cushion, there is an increase against baseline results but as the chart shows, this remains well within the comfort zone. [Fig. 3 & 4]

After ten minutes seating, the Polyair®/GIDL combination continues to further decrease peak pressure, whilst the Integrity® results are only marginally different to the two-minute test results. [Fig. 5 & 6]



MEDICAL CUSHION TYPE:

1. PolyAir® comfort cushion
2. Sumed Integrity® Static High Risk

SUBJECT: Male, 5'6", 82.5kg

DATE OF TEST: April 2021

TESTER: Sumed International (UK) Ltd

*mmHg stands for millimetres of mercury and is used as a pressure measurement

78

42

MmHg* scale

5